




Preparing FACS Aria for NON-Sterile sorting

	<u>STEPS</u>	<u>NOTES</u>
1.	If the computer is off, turn it on, login with: Username - Operator Password - public	
2.	Turn on the FACS Aria	green button on the left side of FACS Aria
3.	Turn on the lasers	White button / mechanical knobs / Software
4.	Upload BD FACSDiva ICON	Diva is a heavy & slow program. Anything within 3 minutes response is normal. Be patient!
5.	Log in with your labs' username and password In the popup window, choose "use CST setting"	For public use, login with Username - public Password - public
6.	Turn on the cooling system ("LAUDA") (beep sound for a couple of seconds is normal)	It takes ~15 minutes to reach default settings of 4°C -8°C . Confirm temperature before you begin to sort.


Checking the fluids in the cart

7.	Check the waste container status. To pass start-up, waste tank should be no more than ¾ full.	For instructions on how to change waste, go to page 9.
8.	Check Sheath status – the stainless steel tank	If needed add ready-made PBS from round 10L white containers the unit holds. (Marked: 'PBS for sorting').
9.	Check Ethanol tank status.	If you find the level to be below the line, add 70% Ethanol. Ready-made 70% Ethanol is at the sink area, in 5 L plastic containers

Fluidic startup

10.	Go to “cytometer” → “fluidic startup”	Instructions will pop up.
11.	Open the sorter cover	
12.	Follow the steps directed by the software with this protocol notes.	After each step you performed press ‘done’
13.		For your safety, <ol style="list-style-type: none"> 1. Validate PBS tank is fully depressurized. 2. First, connect the fluid blue line to the sheath filter. 3. Then connect the transparent air line
14.		Don't be tempted skip this step
15.	 <p>Prepare the nozzle you will use</p>	~5 minutes. At 30%, Airflow enter the sheath tank – at this point Fix air leaks , if any.
16.	After removing the closed nozzle loop and BEFORE you insert the required nozzle: validate that all the slots, sockets, camera area and deflection plates are dry	
17.		
18.	Validate ND filter no 1.	You can find ND filter types in SORTER tool drawer. No 1 is essential for accuracy of RB and AC calibrations

Stream operation

<p>19.</p>	<p>Validate configuration. In the ‘stream window’ you can see configuration set up (picture 1a). If there is a mismatch between the nozzle you inserted and the configuration set-up, go to:</p> <ol style="list-style-type: none"> 1) ‘Cytometer’ → ‘view configuration’ (if a configuration window doesn’t appear, minimize the DIVA, you’ll find the configuration window behind) 2) Choose the correct nozzle configuration 3) Press ‘set configuration’ 4) a message window pops up: Press ‘Ok’. 5) Press ‘Ok’ again at the right side of the window. 6) Close this window at the windows’ X 7) Choose "use CST settings" in the pop up window 	<p>Picture 1 – Stream window</p> 
<p>20.</p>	<p>Turn on the stream and center it</p> <p>***Attention*** beware!!! Don't touch the deflection plates when high voltage warning light is ON.</p>	<p>Click the red X in the “stream window” (picture 1b). It will change to green ✓.</p> <p>For instructions on how to center a stream go to page 9.</p>
<p>21.</p>	<p>Close the deflection plates door</p>	
<p>22.</p>	<p>Close the sorter cover</p>	

23.

Setting the sweet spot

A Table of settings is attached physically on the sorter

Validate that the frequency of the nozzle is according to Table setting values.

Validate that the left value of the gap is according to Table setting values.

Setting drop 1. Drop-1 is the breakoff point of the stream to separated drops (the location of the first separated drop). You can control the breakoff point by changing the amplitude, scroll up or down the amplitude to set the breakoff point. There are two fields for 'drop-1' values: the left shows the last saved value. The right shows current value. Current value will change when the amplitude is changed.


Recommendations for choosing the breaking point:

- Prefer to locate drop-1 breakoff point at the distance range (*i.e.* Drop-1 value) between 200-300 pixels.
- If you have a couple of possible breaking points within that range, then prefer the value, closest to the last saved value. It is not required to reach the exact value, but to find the closest one.
- **Confirm satellite merge.**

Setting the gap - fine tune the amplitude so that the Gap right value match the left value. Once the gap values match (+-2), go to the left field of Drop -1 and type the value you have on the right field of Drop-1 (current).

Activate the 'sweet spot' to keep the stream conditions stable.

Picture 3



Calibrations:

Calibrations can take a few minutes if the system is clean and up to an hour if you find out it needs cleaning. Consider it in your time planning.

Rainbow Calibration

Rainbow beads (RB) are used to check the system performance and that it is sufficiently clean

You should find a 'ready to use' tube in "FACS ONLY solution rack" on the sorter table else you will find additional tubes in the unit refrigerator. The tube is marked RB.

And...If you must,

Prepare new RB tube: transfer 1 drop of RB from the stock vial into ~0.5ml either DW or PBS

- 24) In the Browser window open 'Shared view' (+), **DoubleClick** on the **notebook** named Rainbow or **RB** - If a window pops up, choose "continue".
- 25) Click (+) in the specimen of the current month.
- 26) Insert a new tube and name it with the current day. Activate the new tube (activated tube turns from gray  to  green).
- 27) Go to "cytometer"  "cleaning modes"  **"sample line backflush"**  **"start"**; let it **wash** the sample line for **30 seconds**. Press **"Stop"** and then **"cancel"**.
- 28) **Load** a tube containing FRESH **DW** and let it run for **30 seconds**.

(This step is essential to wash residues from the sample line that destroy the APC signal of the RB)
- 29) **Validate that all the lasers are turned ON**
- 30) **Vortex & Load a tube with RB**
- 31) **Acquire** at a **speed of no more than 200 events/sec and no more than flow rate 2**.
- 32) **Press Record**

33) **Unload** the RB tube and **analyze the result as followed**

- Confirm the gate P1 is set on the correct single bead population
- Confirm that P1 is more than 65%
- Confirm that each one of the histograms has 8 peaks
- Confirm that the CV of the last right peak is no more than 5%

If RB results are satisfactory, proceed to AccuDrop calibration. Do not proceed to AccuDrop calibration if RB calibration is faulty.

Run RB whenever: **1)** Following “start-up” procedure. **2)** If you changed a nozzle. **3)** As part of troubleshooting step.

AccuDrop Calibration

AccuDrop beads (AC) are used to set the optimal drop delay for sorting.

You should find a ‘ready to use’ tube in “FACS ONLY solution rack” on the sorter table else you will find additional tubes in the unit refrigerator. The tube is marked AC.

If you must....Prepare new AC tube: transfer 1 drop of AC from the stock vial into ~0.5ml either DW or PBS

Steps

- 34) Go to ‘Shared view’ **Open ‘AccuDrop /AC’** experiment; by DoubleClicking on the notebook. If a window pop-up, choose “continue”.
- 35) **Press the ‘+’ of folder ‘Global worksheet’**
- 36) **Press ‘+’ of page ‘Global sheet1’ and DoubleClick ‘Sort layout’**
- 37) **Confirm the sort layout conditions:**


Tube device: 4 tubes.

Precision: FINE TUNE.

Target events: Continuous.

The gate “NOT P1” should be set in the left stream position

Notes



- 38) **Press '+' of the specimen and activate the tube.**
- 39) **Click on voltage** (red circle-picture 4) **press 'test sort'** (blue circle – picture 4). Four side streams should appear. The location of the side streams can be adjusted if needed, with their sliders.
- 40) Close side streams by pressing again 'test sort' & close voltage. by pressing again voltage icon.
- 41) **Vortex & Load AC tube**

Assure to have a fit range of events/second:

For 70 micron = 1,500 to 3,000 evt/s

For 85 micron = 800 to 2,000 evt/s


For 100 micron = 600 to 1,500 evt/s

For 130 micron = 400 to 1,200 evt/s

If evt/sec is lower than desired rate, increase the 'flow rate', you can go up as much as flow rate-5 on the other end, if beads are too concentrated **→** dilute them with DW/PBS.

- 42) **Press 'Sort' → Choose 'cancel'.**
- 42.1 **Turn on the voltage;** a side stream appears on the left **→** Confirm that the **Image of the main stream and left-side streams in the 'sort window' are focused.** Do so by turning the AccuDrop laser adjust (Picture 2, silver knob) either clockwise or counterclockwise until the **two** streams are brightest/strongest.
- 42.2 **When streams are focused, press on 'optical filter' icon.**

Picture 4 – Sort window



Side streams sliders – use them to reposition side streams

When you choose 'cancel' beads are sorted and delivered to the waste.

- 42.3 Two squares appear. Each stream should be inside its square.
- 42.4 **Confirm that the right & left sum ~100% altogether.**
- 42.5 Press 'AutoDelay' in the sorting window.

You can also do **manual setting of the drop delay**; in the 'sorting window', gradually increase or decrease the drop delay by clicking its up/down arrows.

Choose the drop delay that yields 95-100% of the events in the left square. For each change you make in the drop delay wait ~2 seconds to let the sorting values stabilize.

- 43) **Unload AC tube**
- 44) properly Close the AC protocol by double clicking its notebook

DIVA 8 randomly deletes gates from experiments that are not closed properly. Opening the next protocol or pressing "Log out" or "quit" is not a proper way to close the experiment.

Run AC whenever: **1)** Following "start-up" and only after RB test passed. **2)** When changing a nozzle. **3)** When 'Drop 1' has changed in more than 20 units. **4)** If you performed "clean the flow-cell" protocol.

Biohazard safety:

Safety level 2 and Human origin cells should Turn on The AMO system (safety regulations to protect people in the sorter area from breathing aerosols)!!! Remember to turn it off when you are done

The system is now ready for **non-sterile** use

Changing waste container:

A full waste container should be replaced with a reserve empty one, placed by the side of the cart. Confirm the empty container has 500ml sodium hypochlorite (EKONOMICA). Leave the full tank near the cart, but if you must. Empty the full waste tank to the sink and add 500 ml of sodium hypochlorite (that can be found under each sink). Place the empty tank by the cart. A full waste container will stop the stream without warning, giving a “clog” message. Keep track of the waste status.

Cleaning a nozzle: put the nozzle in a tube containing 70% ethanol or DW, sonicate for 1 minute and air dry. Sometimes this step has to be repeated several times until the nozzle is clean **Always validate that the nozzle is clean under the microscope**

Centering the stream

*****Attention*** beware!!! Do not touch the deflection plates when high voltage is ON.**

Open the deflection plates’ door. The stream between the two deflection plates should hit the center of a narrow tab (waste drawer). If not, center it by first loosening two adjustment screws and then manually moving the sort box to center position. Re-tighten the screws to lock the position of the box. To handle the screws use an “Allen” type screwdriver (yellow handle), found in the tool drawer of the Aria.

When the stream is centered

