

Micro Writer ML3 Focus Dose Tests

S1805, S1813
small structures
focus-dose calibration

Leonid Tunik
Nanofabrication unit
Chemical research support department

June 6, 2023



מכון
ויצמן
למדע

WEIZMANN
INSTITUTE
OF SCIENCE



מכון
היצמון
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

Dear User,

The following file should help you to find the appropriate parameters to expose your design with MicroWriter.

The entire file contains several charts that describe the effect of major expose parameters: dose (on the x-scale) and focus correction on critical dimensions (CDs) of exposed features on the y-scale. The presented data were collected from the features exposed with the variation of the following parameters:


The applied dose varied between 35 to 75 mJ/cm²

Focus points chosen at autofocus and +/- 0.5 μm

The list of features includes isolated and nested lines and squares, spaces, and circles. Both direct and inverted image structures were exposed and measured. Titles of the charts correspond to the type and dimension of the designed structure.

S1805 photoresist 0.5μm thick on Si was exposed using an x20 lens. 5nm of Cr was deposited for SEM measurements.

S1813 photoresist 1.5μm thick on Si was exposed using an x20 lens. 5nm of Cr was deposited for SEM measurements.





מכון
היצמח
למדע


WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration


Lines

2.5um lines, separated by 2.5um spaces

Positive



Negative









מכון
היצע
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1805 Lines









מכון
היצע
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1813 Lines






מכון
היצמון
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

Spaces

2,3,4,5 μm between the lines of 1.5 μm width









מכון
היצע
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1805 Spaces









מכון
היצע
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1813 Spaces






מכון
היצמון
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

Gaps

2.5 μm and 3 μm between large exposed areas







מכון
הצמח
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1805 Gaps







מכון
היצמח
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1813 Gaps





מכון
היצמ
למדע


WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration


Squares

1.5, 3, 5, 7 μm

Positive



Negative









מכון
היצע
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1805 Squares









מכון
היצע
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1805 Squares cont...









מכון
הצמחה
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1813 Squares









מכון
היצע
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1813 Squares cont...






מכון
היצמון
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

Circles

diameter 1.5 μm






מכון
היצמח
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1805 Circles






מכון
היצע
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1813 Circles





מכון
היצמח
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

Positive vs. Inverted structure

The “golden point” exposure dose could be found to obtain matching final critical dimensions for positive and inverted structures both. Unfortunately, this final critical dimensions value usually will be lower than the designed value. Therefore, to get the desired final critical dimensions, the corresponding value in the design file should be corrected (increased) accordingly.






מכון
היצע
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1805 Positive vs. Inverted “golden point”









מכון
היצמח
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

S1805, S1813 small structures focus-dose calibration

S1813 Positive vs. Inverted “golden point”



* There were difficulties obtaining stable repetitive results on inverted structures (underexposed). Some of the points in the charts are missing.



מכון
ויצמן
למדע

WEIZMANN
INSTITUTE
OF SCIENCE

Micro Writer Focus Dose Tests

Thank you