


# Disposal of non sharp labware containing liquid

Labware (solid waste) must not contain over 2 ml of liquids in each labware, therefore the content should be collected as following:



Empty the liquid into a container with a narrow opening

Place an appropriate sticker and discard at the departmental waste site



Labware with remaining liquid (less than 2ml) should be discarded according to liquid's volatility:

## **Volatile**

Leave the labware open in the chemical hood until full evaporation, then dispose to a regular bin.



## **Non volatile**


Place the labware in a chemical plastic bag. Post a suitable sticker and dispose the bag at the departmental waste site.



**Never** use a biohazard bag to dispose chemical waste!!!



## **Assorted labware that may be disposed in a chemical bag:**



### **YES, YOU CAN**

Eppendorf tubes (liquid content less than 2 ml) may be disposed in the bag with their content. Before disposing, make sure that Eppendorf tubes are tightly capped.



### **NO, YOU CAN'T**

Don't mix in one bag assorted labware together with over 2 ml bottled liquid waste.